

# Mask re-use guidelines

## How to remove your mask to re-use it the same day:

1. **Wash your hands** well with soap and water for at least 20 seconds.

or

Clean your hands well with hand sanitizer with at least 70% alcohol.



2. **Remove your mask:**

- a. Carefully remove your mask by the ties or ear loops. The mask is used, so remove it carefully while leaning forward.

or

- b. Carefully remove your mask by untying lower ties first, then the upper ties. The mask is used, so remove it carefully while leaning forward.



3. **Inspect the mask.** If it is wet, visibly dirty, torn or deformed, request a replacement mask from your manager.

4. **Store the mask** as directed by your manager. This may be in a labeled brown paper bag; on a clean, dry surface; on a paper towel; or in a specific area.



5. **Wash your hands** well with soap and water for at least 20 seconds.

or

Clean your hands well with hand sanitizer with at least 70% alcohol.



# Mask re-use guidelines

## How to put your mask on for re-use:

1. **Wash your hands** well with soap and water for at least 20 seconds.  
or  
Clean your hands well with hand sanitizer with at least 70% alcohol.



2. **Grasp mask:**
  - a. By the ear loops.  
or
  - b. By the upper ties.



3. **Place over your face** and secure the mask behind your ears or around your head (if using ties). Don't touch the front of the mask.



4. **Wash your hands** well with soap and water for at least 20 seconds.  
or  
Clean your hands well with hand sanitizer with at least 70% alcohol.

